MANUFACTURED BY

Glaucoma Cases,  
**Simplified™**

**MicroSurgical Technology, Inc. (MST)**

8415 154th Ave NE

Redmond, WA 98052 United States

Tel: +1.425.861.4002

Toll Free: 1.888.279.3323

www.microsurgical.com

**EUROPEAN UNION AUTHORIZED  
REPRESENTATIVE**

MDSS GmbH

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
30175 Hannover, Germany

**Patent: [www.microsurgical.com/patents](http://www.microsurgical.com/patents)**

P/N 190369-01 Rev. G © 2020 MicroSurgical Technology, Inc.

# TrabEx+™

Goniotomy +I/A Handpiece

	English	Español	Deutsch	Français	Italiano
①	Purge Chamber/ Clear Cap	Cámara de purga/ tapa transparente	Spülkammer/ Durchsichtige Kappe	Chambre de purge/ Bouchon transparent	Camera di spurgo / Cappuccio
②	Raised Indicators For Tip Orientation	Indicadores en relieve para orientación de la punta	Zusätzliche Anzeigen für Ausrichtung der Spitze	Indicateurs en relief pour l'orientation de la pointe	Indicatori per l'orientamento della punta
③	Irrigation Connector	Conector de irrigación	Spüladapter	Raccord d'irrigation	Connettore d'irrigazione
④	Aspiration Connector	Conector de aspiración	Absaugadapter	Raccord d'aspiration	Connettore d'aspirazione
⑤	Smooth Footplate	Placa lisa para pies	Glatte Fußplatte	Pied lisse	Piastra inferiore liscia
⑥	Pointed Tip	Punta Afilada	Scharfe Spitzen	Extrémité pointue	Estremità Appuntita
⑦	Dual Blades	Cuchillas duales	Doppelmesser	Lames doubles	Lame doppie
⑧	Aspiration Port	Puerto de aspiración	Absauganschluss	Orifice d'aspiration	Porta di aspirazione
⑨	Irrigation Port	Puerto de irrigación	Spülanschluss	Orifice d'irrigation	Porta d'irrigazione

## CAUTION

- United States Federal Law restricts this device to sale to or on order of a physician.
- Do not use the "NAME" Handpiece (Handpiece) if you are not trained for intraoperative gonioscopy and ophthalmic angle surgery.
- Read instructions prior to use.
- Contents are sterile provided that the package is not punctured or damaged. Do not use if the package is damaged or open, or if the seal appears damaged.
- Do not use if product is past the expiration date marked on package label.
- Instrument should be removed from the purge chamber and handled in a manner that avoids the tip touching any surfaces prior to its surgical use in the eye. Such tip contact with surfaces may damage or dull the tip.
- Do not use if the instrument contains debris or is bent or otherwise damaged.
- This procedure is performed within the sterile field in an operating room environment.
- The intended users include trained ophthalmic surgeons that perform the surgical procedure, as well as operating room nurses and staff knowledgeable in the setup, support, and use thereof.
- The Handpiece contains stainless steel and some individuals may be allergic or sensitized to metals. Avoid using the Handpiece in such cases.
- The Handpiece is a single use sterile instrument. Do not reuse or re-sterilize. Reused or re-sterilized devices pose risks of infection and/or usage changes to the tip of the instrument that could result in tissue damage.
- After use, the Handpiece is contaminated and can transmit diseases. Properly dispose of the instrument after use in a properly labeled and puncture resistant biohazard sharps container.

### DESCRIPTION AND INTENDED USE:

The Handpiece is a sterile single use ophthalmic knife. It consists of a pencil-like instrument with a graspable white ABS plastic outer shell (~1.5 cm diameter and ~13 cm length), a stainless steel probe (protruding about 2.2 cm from the tip of the white outer shell and covered by a protective purge chamber), irrigation and aspiration tubing extruded from DEHP-free PVC (~15 cm length) terminated with polycarbonate luer connectors. The diameter of the probe insertable into the eye is ~1.1 mm. The tip of the probe features a pointed footplate tip and dual sharpened edges fabricated of stainless steel.

The Handpiece is intended to be used in ophthalmic surgery to manually cut Trabecular Meshwork (TM) in adult and infantile patients

**THERE ARE TWO OPTIONS FOR USE OF THE HANDPIECE (Refer to DIRECTIONS section for proper use of either option):**

- **FIRST OPTION: Use With Ophthalmic I/A System** The Handpiece can be connected to a suitable ophthalmic Irrigation/Aspiration (I/A) system to provide irrigation and aspiration of balanced salt solution (BSS) while the procedure is being performed.
- **SECOND OPTION: Use Without Ophthalmic I/A System** The Handpiece can be used without being connected to a source of I/A.

### GENERAL CONSIDERATIONS:

1. Use of a gonioscopes on the patient's cornea under surgical microscope magnification is necessary for visualizing the TM and the tip of the Handpiece during use.
2. Tilting of the surgical microscope in combination with positioning of the patient's head is essential to attaining a proper view of angle structures.
3. The Handpiece can be used stand-alone or in combination with other surgical procedures.
4. Use 1.8 mm ophthalmic slit knife to properly size corneal incision, which is the minimal size for use with the Handpiece. A 1.8 mm incision is recommended to maintain adequate pressurization and volumetric maintenance of anterior chamber during use of Handpiece with ophthalmic I/A system.
5. The clear corneal incision should be approximately 180° from the TM removal site.
6. The tip of the Handpiece is very delicate and can be damaged from impact or pulling forces. The Handpiece should always be in the purge chamber when it is not in use and should **never** be laid down without the purge chamber in place. Replace the purge chamber when the device is not in use. The purge chamber/clear cap must be removed carefully to prevent damage to the tip.

### PACKAGING AND SHELF-LIFE:

The Handpiece is provided sterile and should be maintained in its original packaging until time of use. Contents are sterile if package is unopened and undamaged. The expiration date is marked on the package label. Do not use if product is past expiration date marked on package label. This product should be stored in a cool, dry environment in its original packaging.

### CONTRAINDICATIONS:

Proper visualization using a Gonioscope is required to manually cut TM with the Handpiece. Do not attempt to use the Handpiece if access to or visualization of the angle structures is not possible, or if corneal clarity is not sufficient.


### WARNINGS:

1. This device should be used by a surgeon who is trained in gonioscopic angle surgery.
2. During ophthalmic procedures, irrigation, viscoelastic, or other chamber maintenance is essential to maintain a safe surgical environment.
3. The tip of the Handpiece must not be manipulated as its structural integrity and positioning could be adversely affected.
4. Do not re-sterilize or reuse any components. Patients can be exposed to sepsis or infection if any single use, sterile component is reused.

### WARRANTY:

MST warrants that its products shall be free from defects in material and workmanship under normal use and service. MST's liability under this warranty is limited to replacing its products. The foregoing warranty excludes and is in lieu of all other warranties whether expressed or implied including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.


MST disclaims any liability for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this product, other than as expressly provided by specific law. MST neither assumes nor authorizes any other person to assume for it any other or additional liability for loss, damage, or expense in connection with this product.

# DIRECTIONS [See Cover Figures]

## FIRST OPTION: Use With Ophthalmic I/A System

Like other anterior chamber eye procedures, when connecting the Handpiece to a suitable ophthalmic I/A system (such as provided on phaco platforms), sufficient irrigation must be provided when aspiration is being used so that the anterior chamber volume is maintained. In general, operating the Handpiece with the BSS bottle or bag elevated at least 40 cm above the patient's surgical eye combined with volumetric (peristaltic mode) aspiration rates not exceeding 5 ml/min are recommended to maintain the anterior chamber. Ophthalmic I/A systems differ – please consult the operation manuals for your I/A system to insure that the settings used provide the adequate levels and balance of irrigation and aspiration.

	English	Español	Deutsch	Français	Italiano
	<b>SYMBOLS USED:</b> Definition of packaging symbols	<b>SÍMBOLOS UTILIZADOS:</b> Definición de los símbolos del envase	<b>VERWENDETE SYMBOLE:</b> Definition der Verpackungssymbole	<b>SYMBOLES UTILISÉS:</b> Définition des symboles d'emballage	<b>SIMBOLI UTILIZZATI:</b> Definizione dei simboli sull'imballaggio
	Manufacturer	Fabricante	Hersteller	Fabricant	Produttore
	Authorized representative in the European Community	Representante autorizado en la Comunidad Europea	Bevollmächtigter Vertreter in der Europäischen Union	Représentant autorisé dans la Communauté Européenne	Rappresentante autorizzato nella Comunità Europea
	"Use by" or expiration date	Fecha de caducidad	Bis Verfallsdatum verwenden	Date d'expiration	"Usare entro" o data di scadenza
	Lot number	Número de lote	Losnummer	Numéro de lot	Numero di lotto
	Model/catalog number	Número de catalogo/modelo	Modell-/ Katalognummer	Numéro de modèle/catalogue	Modello/ Numero di catalogo
	Sterilized using irradiation	Esterilizado mediante radiación	Sterilisiert durch Strahlung	Stérilisé par irradiation	Sterilizzato per irraggiamento
	Sterile	Esterilizado	Steril	Stérile	Sterile
	Do not use if package is damaged	No utilizar si el envase está dañado	Bei beschädigter Verpackung nicht verwenden	Ne pas utiliser si l'emballage est endommagé	Non utilizzare se la confezione è danneggiata
	Do not re-sterilize	No reesterilizar	Nicht wiedersterilisieren	Ne pas restériliser	Non risterilizzare
	Do not re-use	No reutilizar	Nicht wiederverwenden	Ne pas réutiliser	Non riutilizzare
	Consult Instructions for Use	Consultar las instrucciones de uso	Gebrauchsanweisung beachten	Consulter le mode d'emploi	Consultare le istruzioni per l'uso
	Caution	Precaución	Achtung	Mise en garde	Attenzione
	Contents	Contenido	Inhalt	Contenu	Contenuto
	This product does not contain Latex	Este producto no contiene látex.	Dieses Produkt enthält kein Latex	Ce produit ne contient pas de latex	Il prodotto non contiene latex
	Caution: U.S. federal law restricts this device to sale to or on order of a physician	Precaución: la Ley Federal de los EE. UU. restringe la venta de este dispositivo a médicos o a quien tenga la orden de un médico.	Achtung: Nach US-amerikanischem Bundesgesetz darf dieses Gerät nur von einem Arzt oder auf Anordnung eines Arztes verkauft werden	Mise en garde : la loi fédérale des États-Unis limite la vente de cet appareil auprès ou sur prescription d'un médecin	Attenzione: La legislazione degli Stati Uniti limita la vendita di questo dispositivo a o su ordine di un medico

- Inspect all packages prior to use. Do not use if package is open or damaged.
- Remove the Handpiece from its package in the sterile field. If the purge chamber/clear cap (1) is not situated on the Handpiece, the mechanical integrity of the tip may be compromised.
- Connect the irrigation connector (tubing with color stripe) (3) of the Handpiece to the irrigation line connector of the I/A system's fluidics set. Connect the aspiration connector of the Handpiece (4) to the aspiration line connector of the I/A system's fluidics set. Place the Handpiece inclined tip down (with the purge chamber in place).
- Activate the I/A system, and purge the fluidic lines until all bubbles are removed in both the Handpiece and I/A system's irrigation and aspiration lines. During this step, the purge chamber may remain partially filled with air, and/or some BSS may leak from the purge chamber.
- Perform a clear corneal incision using 1.8 mm surgical slit knife. Allowing the anterior chamber to de-pressurize facilitates blood congesting into Schlemm's canal, helping identify the location of the TM.
- Carefully remove the purge chamber/clear cap (1) from the Handpiece so the tip does not touch the purge chamber. The purge chamber is easily removed in the following two-steps:
  - Twist purge chamber first to loosen the cap while holding the handpiece tip inside the purge chamber.
  - Pull the Handpiece tip straight out while avoiding tip contact with the purge chamber.
- The raised indicators on the Handpiece (2) are a guide to the orientation of the Handpiece tip and should be held perpendicular to the 1.8 mm corneal incision. Avoid contacting the corneal endothelium with the Handpiece tip when advancing or retracting the instrument in the anterior chamber.
- Activate irrigation of the I/A system. Insert the tip of the Handpiece through the 1.8 mm corneal incision and into the anterior chamber. The anterior chamber should deepen when the irrigation ports (9) of the Handpiece clear the corneal incision.
- Advance the tip toward the TM opposite of the incision site. Establish gonioscopic view of the TM and proceed with advancing the tip until its footplate contacts the TM.
- Pierce the TM with the pointed tip (6) so it enters Schlemm's Canal while avoiding damage to the back wall.
- Activate aspiration of the I/A system (keeping irrigation also activated). Advance the tip of the Handpiece along Schlemm's Canal in the direction you intend to cut TM. Keep the smooth footplate (5) moving along Schlemm's canal while only exerting minimal pressure against the back wall of Schlemm's canal.
 

**CAUTION:** Maintain both dual blades (7) perpendicular to the TM during TM excision or only one of the dual blades may cut the TM.
- As TM is excised, blood reflux from collector channels may appear. The irrigation(9)/aspiration(8) functions of the Handpiece improve visualization of the treatment area.
- The aspiration function of the handpiece removes the excised TM strip if the end of the strip is amputated. If the end of TM strip is not amputated, then a micro-forceps can be used to manually grasp and remove the TM strip.
- Gently remove the tip of the Handpiece from the corneal incision following the same orientation as for insertion.
- At the end of the procedure, it is recommended to thoroughly remove any viscoelastic that may have been used, using standard I/A practice to avoid post-operative pressure spikes and place a suture.

## SECOND OPTION: Use Without Ophthalmic I/A System

- Inspect all packages prior to use. Do not use if package is open or damaged.
- Remove the Handpiece from its package in the sterile field. If the purge chamber is not situated on the Handpiece, the mechanical integrity of the tip may be compromised.
- It is important to remove air from the Handpiece and its I/A pathways to avoid air bubbles entering the anterior chamber. Connect the Handpiece irrigation connector (tubing with color stripe) (3) to a syringe containing BSS. Place the Handpiece inclined tip down with attached purge chamber/clear cap (1) and hold in place. Infuse BSS until the BSS exits the aspiration connector (4). Remove the syringe and connect the irrigation connector (3) and the aspiration connector (4) of the Handpiece together. During this step, the purge chamber may remain partially filled with air, and/or some BSS may leak from the purge chamber.
- Perform a clear corneal incision using 1.8 mm surgical slit knife. Allowing the anterior chamber to de-pressurize facilitates blood congesting into Schlemm's canal, helping identify the location of the TM.
- Inject a viscoelastic into the anterior chamber to inflate and maintain its volume during the procedure.
 

**CAUTION:** The volume of anterior chamber must be re-established and maintained per this step for continuing the procedure.
- Carefully remove the purge chamber/clear cap (1) from the Handpiece so the tip does not touch the purge chamber. The purge chamber is easily removed in the following two-steps:
  - Twist purge chamber first to loosen the cap while holding the handpiece tip inside the purge chamber.
  - Pull the Handpiece tip straight out while avoiding tip contact with the purge chamber
- The raised indicators on the Handpiece (2) are a guide to the orientation of the Handpiece tip and should be held perpendicular to the 1.8 mm corneal incision. Avoid contacting the corneal endothelium with the Handpiece tip when advancing or retracting the instrument in the anterior chamber.
- Insert the tip of the Handpiece through the 1.8 mm corneal incision and into the anterior chamber.
- Advance the tip toward the TM opposite of the incision site. Establish gonioscopic view of the TM and proceed with advancing the tip until its footplate contacts the TM.
- Pierce the TM with the pointed tip (6) so it enters Schlemm's Canal while avoiding damage to the back wall.
- Advance the tip of the Handpiece along Schlemm's Canal in the direction you intend to cut TM. Keep the smooth footplate (5) moving along Schlemm's canal while only exerting minimal pressure against the back wall of Schlemm's canal.
 

**CAUTION:** Maintain both dual blades (7) perpendicular to the TM during TM excision or only one of the dual blades may cut the TM.
- As TM is excised, blood reflux from collector channels may appear. Use of viscoelastic and/or BSS are methods for improving visualization.
- The TM strip should not be left in the anterior chamber. A micro-forceps can be used to manually grasp and remove the TM strip.
- Gently remove the tip of the Handpiece from the corneal incision following the same orientation as for insertion.
- At the end of the procedure, it is recommended to thoroughly remove any viscoelastic that may have been used, using standard I/A practice to avoid post-operative pressure spikes and place a suture.

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